

## The NCAR Trace Organic Gas Analyzer with Time of Flight mass spectrometer (TOGA-TOF)

The TOGA-TOF is an in-situ online gas chromatograph mass spectrometer (GC-MS) that provides near-continuous atmospheric mixing ratios of an extensive and growing list of volatile organic compounds (VOCs) in the C<sub>1</sub>-C<sub>10</sub> molecular structure range. The list of species that can be quantified using TOGA-TOF includes alkanes, alkenes, aromatics, halogenated VOCs, nitrates, nitriles, sulfides, alcohols, ketones, aldehydes, esters, and ethers. One hundred or more unique trace gases can be measured, with sufficient sensitivity to measure trace species in the remote background atmosphere and dynamic range to measure in highly polluted regions. Table 1 is a list of species that have been identified and/or quantified using the TOGA-TOF, as well as typical limits of detection where known.

The TOGA-TOF uses a cryogenic preconcentrator, consisting of a custom-built liquid nitrogen (LN<sub>2</sub>) dewar and a system of traps, for water removal, sample enrichment, and cryofocusing of trace constituents. A high-resolution electron impact time-of-flight mass spectrometer (HR-EI-TOF-MS; Tofwerk, Switzerland) is used to identify and quantify individual compounds, using Igor-based software packages developed by Tofwerk and Aerodyne Research Inc. System calibrations and blanks are performed using a catalytic-clean air generator/dynamic dilution system with accurate (±1%) and precise (±1%) calibration gas delivery. The system operates continuously, allowing for frequent calibrations and zeros during flight. The TOGA-TOF inlet is a recently-designed constant mass flow design to mitigate the intrusion of aerosols and allow for calibrations and system blanks through the inlet and associated tubing. The TOGA-TOF is contained in an elongated HIAPER rack standard DC-8 rack, weighs less than 200 kg and consumes ~1 kW of power. More information on the TOGA and TOGA-TOF can be found here: <https://www2.acom.ucar.edu/voc-measurements/measurement-instrumentation>.

**Table 1.** Chemical species identified using the TOGA-TOF and limits of detection (LOD) where known.

NMHCs	Formula	LOD, ppt	OVOCs	Formula	LOD, ppt	OVOCs	Formula	LOD, ppt	Halogenated VOCs	Formula	LOD, ppt
<b>Alkanes</b>			<b>Aldehydes</b>			<b>Ethers/Furans (cont'd)</b>			<b>Halogenated VOCs</b>		
Propane	C <sub>3</sub> H <sub>8</sub>	5	Formaldehyde	CH <sub>2</sub> O	20	2-Methylfuran	C <sub>5</sub> H <sub>8</sub> O	0.5	Methyl chloride	CH <sub>3</sub> Cl	1
Isobutane	C <sub>4</sub> H <sub>10</sub>	1	Acetaldehyde	C <sub>2</sub> H <sub>4</sub> O	5	3-Methylfuran	C <sub>5</sub> H <sub>8</sub> O	0.5	Dichloromethane	CH <sub>2</sub> Cl <sub>2</sub>	1
<i>n</i> -Butane	C <sub>4</sub> H <sub>10</sub>	1	Propanal	C <sub>3</sub> H <sub>6</sub> O	2	2,3-Dimethylfuran	C <sub>6</sub> H <sub>8</sub> O	TBD	Chloroform	CHCl <sub>3</sub>	1
Isopentane	C <sub>5</sub> H <sub>12</sub>	1	Isobutanal	C <sub>4</sub> H <sub>8</sub> O	TBD	2,4-Dimethylfuran	C <sub>6</sub> H <sub>8</sub> O	TBD	Tetrachloromethane	CCl <sub>4</sub>	1
<i>n</i> -Pentane	C <sub>5</sub> H <sub>12</sub>	1	Butanal	C <sub>4</sub> H <sub>8</sub> O	1	2,5-Dimethylfuran	C <sub>6</sub> H <sub>8</sub> O	TBD	Tetrachloroethane	C <sub>2</sub> Cl <sub>4</sub>	0.1
2-Methylpentane	C <sub>6</sub> H <sub>14</sub>	0.5	Acrolein	C <sub>3</sub> H <sub>4</sub> O	1	2-Ethylfuran	C <sub>6</sub> H <sub>8</sub> O	TBD	1,2-Dichloroethane	C <sub>2</sub> H <sub>4</sub> Cl <sub>2</sub>	0.5
3-Methylpentane	C <sub>6</sub> H <sub>14</sub>	0.5	Methacrolein	C <sub>4</sub> H <sub>6</sub> O	1	3-Ethylfuran	C <sub>6</sub> H <sub>8</sub> O	TBD	Methyl chloroform	CH <sub>2</sub> Cl <sub>3</sub>	0.5
<i>n</i> -Hexane	C <sub>6</sub> H <sub>14</sub>	0.5	2-Butenal	C <sub>4</sub> H <sub>6</sub> O	2	2-Vinylfuran	C <sub>6</sub> H <sub>8</sub> O	TBD	Methyl bromide	CH <sub>3</sub> Br	1
<i>n</i> -Heptane	C <sub>7</sub> H <sub>16</sub>	1	Furfural	C <sub>5</sub> H <sub>4</sub> O <sub>2</sub>	TBD	3-Vinylfuran	C <sub>6</sub> H <sub>8</sub> O	TBD	Dibromomethane	CH <sub>2</sub> Br <sub>2</sub>	0.03
2,2,4-Trimethylpentane	C <sub>8</sub> H <sub>18</sub>	0.5	3-Furaldehyde	C <sub>5</sub> H <sub>4</sub> O <sub>2</sub>	TBD	1,3-Dioxolane	C <sub>3</sub> H <sub>6</sub> O <sub>2</sub>	TBD	Bromoform	CHBr <sub>3</sub>	0.1
<i>n</i> -Octane	C <sub>8</sub> H <sub>18</sub>	0.5	<b>Ketones</b>			<b>Nitrogen-containing VOCs</b>			Methyl iodide	CH <sub>3</sub> I	0.03
<b>Alkenes</b>			Acetone	C <sub>3</sub> H <sub>6</sub> O	5	<b>Formula</b>			Diodomethane	CH <sub>2</sub> I <sub>2</sub>	0.05
Propene	C <sub>3</sub> H <sub>6</sub>	5	MEK	C <sub>4</sub> H <sub>8</sub> O	0.5	<b>LOD, ppt</b>			Ethyl iodide	C <sub>2</sub> H <sub>5</sub> I	0.5
1-Butene/Isobutene	C <sub>4</sub> H <sub>8</sub>	1	MVK	C <sub>4</sub> H <sub>6</sub> O	0.5	<b>Nitriles</b>			Chloriodomethane	CH <sub>2</sub> ClI	0.05
<i>cis</i> -2-Butene	C <sub>4</sub> H <sub>8</sub>	1	2,3-Butanedione	C <sub>4</sub> H <sub>6</sub> O <sub>2</sub>	TBD	Hydrogen cyanide	HCN	5	Bromodichloromethane	CHBr <sub>2</sub> Cl	0.05
<i>trans</i> -2-Butene	C <sub>4</sub> H <sub>8</sub>	1	<b>Alcohols</b>			Chlorine cyanide	ClCN	TBD	Dibromochloromethane	CHBr <sub>2</sub> Cl	0.03
Isoprene	C <sub>5</sub> H <sub>8</sub>	1	Methanol	CH <sub>3</sub> O	5	Acetonitrile	C <sub>2</sub> H <sub>3</sub> N	1	CFC-11	CCl <sub>3</sub> F	5
α-Pinene	C <sub>10</sub> H <sub>16</sub>	1	Ethanol	C <sub>2</sub> H <sub>5</sub> O	2	Propanenitrile	C <sub>3</sub> H <sub>5</sub> N	1	CFC-12	CCl <sub>2</sub> F <sub>2</sub>	1
β-Pinene/Myrcene	C <sub>10</sub> H <sub>16</sub>	1	2-Propanol	C <sub>3</sub> H <sub>7</sub> O	4	Acrylonitrile	C <sub>3</sub> H <sub>3</sub> N	1	CFC-113	C <sub>2</sub> Cl <sub>3</sub> F <sub>3</sub>	1
Camphene	C <sub>10</sub> H <sub>16</sub>	1	Ethanol	C <sub>2</sub> H <sub>5</sub> O	TBD	Methylacrylonitrile	C <sub>4</sub> H <sub>5</sub> N	2	CFC-114	C <sub>2</sub> Cl <sub>2</sub> F <sub>4</sub>	1
Limonene/3-Carene	C <sub>10</sub> H <sub>16</sub>	1	MBO (2-Methyl-3-Buten-2-ol)	C <sub>5</sub> H <sub>10</sub> O	1	<b>Nitrates</b>			HCFC-22	CHClF <sub>2</sub>	1
Tricyclene	C <sub>10</sub> H <sub>16</sub>	1	<b>Esters</b>			Methyl nitrate	CH <sub>3</sub> NO <sub>3</sub>	TBD	HCFC-141b	C <sub>2</sub> H <sub>3</sub> Cl <sub>2</sub> F	1
<b>Aromatics</b>			Methyl formate	C <sub>2</sub> H <sub>4</sub> O <sub>2</sub>	TBD	Ethyl nitrate	C <sub>2</sub> H <sub>5</sub> NO <sub>3</sub>	TBD	HCFC-142b	C <sub>2</sub> H <sub>3</sub> ClF <sub>2</sub>	TBD
Benzene	C <sub>6</sub> H <sub>6</sub>	0.3	Methyl acetate	C <sub>3</sub> H <sub>6</sub> O <sub>2</sub>	TBD	Isopropyl nitrate	C <sub>3</sub> H <sub>7</sub> NO <sub>3</sub>	2	HFC-134a	C <sub>2</sub> H <sub>2</sub> F <sub>4</sub>	1
Toluene	C <sub>7</sub> H <sub>8</sub>	0.3	Methyl propionate	C <sub>4</sub> H <sub>8</sub> O <sub>2</sub>	TBD	<i>n</i> -Propyl nitrate	C <sub>3</sub> H <sub>7</sub> NO <sub>3</sub>	2	<b>Sulfur-containing VOCs</b>		
Ethylbenzene	C <sub>8</sub> H <sub>10</sub>	0.2	Ethyl acetate	C <sub>4</sub> H <sub>8</sub> O <sub>2</sub>	TBD	<i>t</i> -Butyl nitrate	C <sub>4</sub> H <sub>9</sub> NO <sub>3</sub>	2	Carbonyl sulfide	OCS	1
<i>p</i> -/m-Xylene	C <sub>8</sub> H <sub>10</sub>	0.2	<b>Ethers/Furans</b>			2-Butyl+isobutyl nitrate	C <sub>4</sub> H <sub>9</sub> NO <sub>3</sub>	2	Carbon disulfide	CS <sub>2</sub>	0.2
<i>o</i> -Xylene	C <sub>8</sub> H <sub>10</sub>	0.2	MTBE (Methyl <i>t</i> -butyl ether)	C <sub>5</sub> H <sub>12</sub> O	0.3	<i>n</i> -Butyl nitrate	C <sub>4</sub> H <sub>9</sub> NO <sub>3</sub>	2	Methanethiol	CH <sub>3</sub> S	TBD
Styrene	C <sub>8</sub> H <sub>8</sub>	0.1	Furan	C <sub>4</sub> H <sub>6</sub> O	1	<b>Other</b>			DMS (Dimethyl sulfide)	C <sub>2</sub> H <sub>6</sub> S	1
Ethylbenzene	C <sub>8</sub> H <sub>8</sub>	TBD	THF (Tetrahydrofuran)	C <sub>4</sub> H <sub>8</sub> O	TBD	Pyrrrole	C <sub>4</sub> H <sub>5</sub> N	TBD	Carbon suboxide	C <sub>3</sub> O <sub>2</sub>	TBD
						Nitromethane (CH <sub>3</sub> NO <sub>2</sub> )	CH <sub>3</sub> NO <sub>2</sub>	TBD			

**Measurements:**

VOCs

**Aircraft:**

Gulfstream V – NSF, C-130 – NSF, DC-8 - AFRC

**Measurement Characteristics:**

Overall estimate of uncertainty: VOC dependent, typically 20%.

**Precision:** 3% or less.

**Response time:** typically 35-second samples on a 2-minute cycle, with the option for a shorter, targeted sample.

**Instrument Type:**

Gas chromatography, Time-of-Flight Mass Spectrometry (in situ)

**Point(s) of Contact:**

Eric Apel (POC; PI)